

NEWLIFE® INTENSITY®

Stationary Oxygen Concentrator

HIGH FLOW AND HIGH PRESSURE COMBINED

Combining high flow and high output pressure makes the NewLife Intensity the premier high flow stationary concentrator for at home and facility use.



HIGH-PRESSURE APPLICATIONS

The Intensity has an output pressure of 20 psi (137.9 KPa), more than double any other high flow concentrator.



HIGH FLOW DELIVERY

Offering up to 10 LPM, this time-proven concentrator offers higher flow options to serve the patients in hospitals, long-term care facilities, clinics, waiting rooms, or day and night time settings.



DUAL FLOW OPTION

The unique Dual Flow option features two separate flow meters that can provide a total of 10 LPM — allowing two patients to get oxygen from one machine simultaneously.

CAIRE Contact Information:

Providers Toll Free

1.800.482.2473

Patients Toll Free

1.877.704.0878

www.caireinc.com



SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© Copyright 2020 CAIRE Inc. All Rights Reserved.

NEWLIFE® INTENSITY®

Stationary Oxygen Concentrator

Includes

- Oxygen Concentrator
- Cannula Outlet Fitting

Product Features

- High output pressure
- Time-proven NewLife ruggedness and reliability
- Oxygen monitor
- Dual flow option available
- Easy to service and maintain
- Dual function hour meter

Specifications

FLOW RATE
2-10 LPM

O₂ CONCENTRATION
90%+5.5%/-3%

MAXIMUM TUBING
200 ft (61 m)

DIMENSIONS
27.5 x 16.5 x 14.5 in (69.9 x 41.9 x 36.8 cm)

OPERATING ENVIRONMENT
41° to 104°F (5° to 40°C),
15-90% humidity

MAINTENANCE SCHEDULE
Felt Filter – 1 Year Replacement
Intake Filter – Clean Weekly

WEIGHT
58 lbs (26.3 kg)

STORAGE ENVIRONMENT
-13 to 158°F (-25 to 70°C)
0-90% humidity

WARRANTY
2 Years

SOUND LEVEL
58 dB(A)



Available Accessories



HUMIDIFIER BOTTLE

Ask your doctor or therapist if humidification is right for you.



NASAL CANNULAS

Available in 7 ft (2.1 m), 25 ft (7.6 m).
Talk to our Customer Service team to learn more.



OXYGEN TUBING

Available in 7 ft (2.1 m), 25 ft (7.6 m), and 50 ft (15.2 m).

CAIRE Contact Information:

Providers Toll Free

1.800.482.2473

Patients Toll Free

1.877.704.0878

www.caireinc.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© Copyright 2020 CAIRE Inc. All Rights Reserved.

ML-CONC0084 G